
High Throughput Assays for the Characterization of Monomeric GTPases of the Ras Superfamily


Michael Lassle¹, Marie-Elaine Caruso², Sarah Jenna², Simon Beaulne³, Eun Hye Lee², Anne Bergeron⁴, Cédric Chauve³, Philippe Roby¹, Jean-François Rual⁵, David E Hill⁵, Marc Vidal⁵, Roger Bossé¹ and Eric Chevet². ¹PerkinElmer BioSignal (Canada), ²McGill University (Canada), ³U.Q.A.M. (Canada), PerkinElmer' 🐯 McGill precisely. ⁴University of Montreal (Canada), ⁵Dana-Farber Cancer Institute (USA)

Abstract


We have developed an innovative analytical platform to perform systematic comparison of protein families based on their biochemical properties. Our approach was validated on the Rho subfamily of GTPases. We used two high throughput technologies (AlphaScreen[™] and FlashPlate®), to measure nucleotide binding capacity, exchange, and hydrolysis activities of small monomeric GTPases. These two technologies benefit from high sensitivity while allowing for homogenous and high throughput assays. To analyze and integrate the data obtained, we have developed an algorithm that allows the classification of GTPases according to their enzymatic activities. Integration and hierarchical clustering of these results revealed unexpected features of the small Rho GTPases when compared with primary sequence-based trees. We therefore suggest a novel phylobiochemical classification of the Ras superfamily of GTPases



A) SDS-PAGE of RhoGTPases GST-fusion proteins. Cloning was performed using Gateway from Invitrogen. B) A Rho GTPases phylogenetic tree was produced by multiple sequence alignments a with ClustaW followed by DNAsis. All GTPases except RhoBH and RacI M are found within the




[35S]GTP/S saturation binding assay. FlashPlate® experimental design for [35S]GTPyS binding assay, Nucleotide-free GST-GTPases were added to a GSH FlashPlate and incubated with different concentrations of [35S]GTPvS. Radioligand interacting with the GST-GTPase activates the FlashPlate scintillation and generate a luminescence




the [35S]GTPγS and GST-CDC-42 or GST alone. The apparent Kd (Kd_{ann}) value corresponds to (a) Control and Control of Con GTPase (mRAC-1N17). All the Kdapp values represent the average of at least three indepen iments performed in trinli

5 **GTP** competition activity infin-GTP/S competition assay. AlphaScreen experimental design for the hiotin-GTP/S


oun-of 1/5 competition assay. Application experimental design to the order of 1/5 impetition assay. Nucleotide-free GST-GTPases were incubated with a constant concentration otin-GTP/S and increasing concentrations of unlabeled competitors. Streptavidin Donor beads cognize biotin-GTP/S and the anti-GST Accentor heads bind to GST-GTPases. When biotinecognize boom-G1P/S, and the ann-GS1 Acceptor beads bind to GS1-G1Pases. When boom-TP/S interacts with the protein, the two beads are in close proximity and a signal is detectable pon laser excitation of the Donor bead. Presence of unlabeled nucleotides (GDP, GTP) compete he binding of the biotin-GTP/S tracer to the GST-GTPase producing a signal decrease.




a) biotin-GTPyS competition curves obtained with either ATP or GTP on CDC-42 or GST. The curves and the result of a single experiment done in triplicate. The IC₃₀ value corresponds to the competitor concentration required for the loss of 50% of the maximum signal (obtained in the absence of competitor) B) comparison of the IC₃₀ values obtained for six C elegans GTPases including a GTP binding-deficient TPase (mRAC-1N17) using the biotin-GTP₁S competition assay. The IC₅₀ values for the competition ssay represent the average of at least three independent experiments performed in triplicate.



frolysis activity. GST-GTPases were added to a GSH FlashPlate and were loaded with eithe 500 nCi of [35SIGTPyS or 500 nCi of [33P]GTP. When the GTPase hydrolyzes the y phosphate 13 PJGTP, the radioactivity is liberated in the reaction mixture, and β particles are quenched by the aqueous solution leading to the inactivation of the scintillation reaction. A decrease of signal is then measured. [35S]GTP/S, a non-hydrolyzable analog of the GTP, is used as a control to measure the nucleotide dissociation from the GST fusion protein.





based on GTPase biochemical activities are shown (A-D) and according to the integration of these four activities (E). These trees were then compared with the amino acid sequence- based tree (F). The bar graph (G) shows the distance, calculated using Euclidean matrices, between the SB tree and different one, two, three, or four activities (A, B, C, and D). The



omology, all GTPases have distinctive properties based on their biochemical activities

Results obtained allowed to propose a new classification of small G-proteins based on their nzymatic properties which provides more pertinent information directly linked to protein

Using the methods described herein, we are currently characterizing all the 61 C.elegans GTPases at once. All proteins will be compared and classified according to their inf biochemical activities. A new cluster tree will be built accordingly.

33PlGTP hydrolys ctivity assay. FlashPlat 8 **GTP hydrolysis activity**

lo ô 🗐 100 4 40° 100° 40° 40° 40° 10° A) the half-life value is the time needed for a GTPase to hydrolyze 50% of the initially bound [33P]GTP with



biotin-GDP/GTP exchange activity assay. AlphaScreen experimental design for the measur of the GTPase exchange activity. GST-GTPases were loaded with biotin-GDP and incubated with different concentrations of GTP. Streptavidin Donor beads capture biotin-GDP while the anti-GST Acceptor beads bind to GST-GTPases. When biotin-GDP interacts with the GTPases the two beads are brought in close proximity and a signal is produced upon laser excitation of the Donor beads. When biotin-GDP is exchanged to GTP, the distance between the beads are no longer in proximity and a loss of signal is measured.

PerkinElmer Life and Analytical Sciences, 45 William Street, Wellesley, MA USA (800) 762-4000 (+1) 203 925-4602 www.perkinelmer.com