

Europium-labeled Ligands for Receptor Binding Studies



Katja Sippola, Maija-Liisa Mäkinen, Sanna Rönmark, Joni Helenius, Riitta Heilimö, Anu Koivikko, Pertti Hurskainen and Christel Gripenberg-Lerche
PerkinElmer Life and Analytical Sciences, Wallac Oy

1

Introduction

Time-resolved fluorescence enhancement technique DELFIA® enables development of highly sensitive assays for screening. We have developed a family of Europium-labeled peptides and proteins designed for ligand receptor binding assays that can be easily automated and optimized either for 96 or 384 well format. These Eu-labeled ligands provide an excellent non-radioactive alternative that is both stable and sensitive. The new tachycinin members of the DELFIA® ligand family are Substance P and Neurokinin A. These peptide ligands have many physiological roles, e.g. they stimulate smooth muscle contraction and glandular secretion and are involved in immune responses and neurotransmission.


2

Methods


The DELFIA® ligand receptor binding assay is based on dissociation-enhanced time-resolved fluorescence. DELFIA® Eu-labeled ligand and receptor membrane prepare are incubated on a filter plate (PALL AcroWell or AcroPrep) after which unbound labeled ligand is removed by filtration. Eu is dissociated from the bound ligand by using DELFIA® Enhancement Solution. Dissociated Eu creates highly fluorescent complexes, which are measured in a multilabel counter with TRF option, e.g. EnVision™.

3

DELFIA® Neurokinin A assays on 96 well filtration plates



The saturation experiment was performed on a 96 well filtration plate with increasing amounts of DELFIA Eu-labeled neurokinin A in the presence of 1 µg of human NK-2 receptor (B_{max} 2.8 pmol/mg protein) per well. Non-specific binding was determined in the presence of 250 nmol/L unlabeled human neurokinin A. A K_d value of 1.7 nmol/L was obtained.




The competition between DELFIA Eu-labeled neurokinin A and unlabeled human neurokinin A. The displacement curve was performed with 1.7 nmol/L of DELFIA Eu-labeled neurokinin A and increasing amounts of unlabeled neurokinin A in the presence of 1 µg of human NK-2 receptor (B_{max} 2.8 pmol/mg protein) per well. A K_i value of 0.4 nmol/L was obtained.

S/B 12
 z' 0,73
 K_d 1,7 nM
 K_i 0,4 nM


Results are comparable to radioligand binding assays

4

DELFIA® Substance P assays on 96 well filtration plates



The saturation experiment was performed on a 96 well filtration plate with increasing amounts of DELFIA Eu-labeled substance P in the presence of 12 µg of human endogenous NK-1 receptor (B_{max} 0.25 pmol/mg protein) per well. Non-specific binding was determined in the presence of 500 nmol/L unlabeled human substance P. A K_d value of 0.61 nmol/L was obtained.



The competition between the DELFIA Eu-labeled substance P and unlabeled human substance P is shown in Figure 2. The displacement curve was performed with 0.7 nmol/L of DELFIA Eu-labeled substance P and increasing amounts of unlabeled substance P in the presence of 12µg of human endogenous NK-1 receptor (B_{max} 0.25pmol/mg protein) per well. A K_i value of 0.34 nmol/L was obtained.

S/B 9
 z' 0,67
 K_d 0,61 nM
 K_i 0,34 nM

Results are comparable to radioligand binding assays

5

Assays on 384 well filtration plates

- 5 µl unlabeled ligand
- 5 µl Eu-ligand
- 10 µl receptor
- incubation 90min (15 s shake)
- filtration wash 3x300 µl (DELFIA® L*R wash solution)
- 20 µl Enhancement Solution
- incubation 15 min (shaking)
- TRF measurement

Kd and Ki values on 384 AcroPrep plates

	Kd	Ki	S/B
Neurokinin A	0,94	0,45	9
Substance P	0,43	0,31	9

Results are comparable to results on 96 well plates

6

Comparison of manual and automated assay formats

Manual protocol

- 25 µl unlabeled ligand
- 25 µl Eu-ligand
- 50 µl receptor
- incubation 90min (15 s shake)
- filtration wash 3x300 µl
- 200 µl Enhancement Solution
- incubation 15 min (shaking)
- TRF measurement

Automation protocol

- 2 µl unlabeled ligand in DMSO
- 50 µl Eu-ligand
- 50 µl receptor
- incubation 2 h - 5 h
- filtration wash 3x300 µl
- 200 µl Enhancement Solution
- incubation 90 min - 4 h
- TRF measurement

	automated assay		manual assay	
	S/B	z'	S/B	z'
Neurokinin A	6,3	0,73	7,2	0,69
Substance P	6,9	0,76	6,7	0,64

Assays can be easily automated - increased through-put with less variation

7

DELFIA® Eu-ligand products available

Motilin	IL-8
Galanin	IL-2
EGF	IL-5
Neurotensin	TNF α
Neurokinin A	Bombesin
Substance P	NDP- α MSH

L*R binding buffer concentrate (10 x)
 L*R wash concentrate (25 x)

8

Conclusions

The DELFIA® ligand receptor binding assay characteristics:

- High sensitivity** - works with low expression level receptors, endogenous receptors
- Non-radioactive** - no radioactive waste, long shelf life
- Several assay formats** - manual and automated, 96- and 384-well format