

AlphaLISA JMJD2A Histone H3-Lysine 36 Demethylase Assay

AlphaLISA®

AlphaLISA #8

Authors

Claire Normand
Mathieu Arcand
Julie Blouin
Mireille Caron
Anne Labonté
Lucille Beaudet
Jaime Padrós

PerkinElmer, Inc.
Montreal, QC
Canada, H3J 1R4

This AlphaLISA immunodetection assay measures the demethylation of a biotinylated Histone H3 (21-44) peptide tri-methylated at lysine 36.

Anti-di-methyl-Histone H3 Lysine 36 (H3K36me2) AlphaLISA® Acceptor Beads

- AL123C: 250 µg, 500 assay points*
- AL123M: 5 mg, 10,000 assay points*
- AL123R: 25 mg, 50,000 assay points*

*0.5 µg/assay point

Peptidic Substrate Sequence:

ATKAARKSAPATGGVK(me3)KPHRYRP-GG-K(Biotin)-OH

AlphaLISA Assays

AlphaLISA technology is a powerful and versatile platform that offers highly sensitive, no-wash immunoassays using Alpha Donor and AlphaLISA Acceptor beads. In this technical note, we present the optimization of an epigenetic enzymatic assay using a biotinylated Histone H3-derived peptide as substrate. Detection of the modified substrate was performed by the addition of Streptavidin (SA) Alpha Donor beads and AlphaLISA Acceptor beads conjugated to an antibody (Ab) directed against the epigenetic mark of interest. Upon laser irradiation of the beads-target complexes at 680 nm, short-lived singlet oxygen molecules produced by the Donor beads can reach the Acceptor beads in proximity to generate an amplified chemiluminescent signal at 615 nm. The intensity of light emission is proportional to the level of biotinylated substrate modification.




Figure 1. Schematic representation of AlphaLISA detection of a modified histone peptide.


Development of a JMJD2A Histone H3-Lysine 36 Demethylase Assay:

Reagents needed for the assay:

Anti-di-methyl-Histone H3 Lysine 36 (H3K36me2) AlphaLISA Acceptor beads	PerkinElmer # AL123
Alpha Streptavidin Donor beads	PerkinElmer # 6760002
Histone H3 (21-44), H3K36(me3) peptide, biotinylated	AnaSpec # 64441
AlphaLISA 5X Epigenetics Buffer 1 Kit	PerkinElmer # AL008
JMJD2A (human), recombinant	BPS BioScience # 50103
White opaque OptiPlate™-384	PerkinElmer # 6007299
TopSeal™-A films	PerkinElmer # 6005185
α-Ketoglutaric acid potassium salt (2OG)	Sigma # K2000
(+) Sodium L-ascorbate	Sigma # 11140
Ammonium iron(II) sulfate hexahydrate (Fe(II))	Sigma # 215406
2,4-Pyridinedicarboxylic acid (2,4-PDCA)	Sigma # P63395


Assay Buffer: 50 mM Hepes pH 7.5, 0.01% Tween-20 and 0.1 % BSA.

Experiment 1: Enzyme Titration and Time-Course



Enzymatic progress curves were performed by incubating JMJD2A at concentrations ranging from 0.125 to 1 nM with 100 nM biotinylated Histone H3K36me3 peptide substrate plus 100 μM 2OG, 5 μM Fe(II) and 100 μM ascorbate. Acceptor beads were added at the indicated times. Donor beads were added 60 min later and signal was read after 30 min. A 30 min reaction time using 0.5 nM enzyme was selected for all subsequent experiments

Experiment 2: 2OG Titration



Serial dilutions of 2OG ranging from 10 nM to 100 μM were added to 0.5 nM JMJD2A and 100 nM biotinylated Histone H3K36me3 peptide substrate plus 5 μM Fe(II) and 100 μM ascorbate. A 2 μM 2OG concentration was selected for subsequent experiments.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com


Standard Protocol

- Dilute JMJD2A enzyme, 2OG, Fe(II), ascorbate, 2,4-PDCA (inhibitor) and biotinylated peptide substrate in Assay Buffer just before use.
- Add to the wells of a white OptiPlate-384:
 - 5 μL of inhibitor (2X) or Assay Buffer
 - 2.5 μL of enzyme (4X)
 - 2.5 μL of biotinylated Histone H3K36me3 peptide/2OG/Fe(II)/ascorbate mix (4X).

For 2OG titration, add 2OG dilutions independently of substrate.
- Cover the plate with TopSeal-A film and incubate at room temperature (RT).
- Prepare 1X Epigenetics Buffer 1 as recommended in the buffer technical data sheet.
- Prepare a 5X Acceptor beads solution at 100 μg/mL in 1X Epigenetics Buffer 1 (final concentration of 20 μg/mL in 25 μL total assay volume).
 - 5 μL of Acceptor beads


Addition of Acceptor beads prepared in 1X Epigenetics Buffer 1 stops the enzymatic reaction.
- Cover with TopSeal-A film and incubate for 60 min at RT.
- Prepare a 2.5X Streptavidin Donor beads solution at 50 μg/mL in 1X Epigenetics Buffer 1 (final concentration of 20 μg/mL in 25 μL total assay volume) in subdued light.
 - 10 μL of Streptavidin Donor beads
- Cover with TopSeal-A film and incubate in subdued light for 30 min at RT.
- Read signal in Alpha mode with the EnVision® or EnSpire® reader.

Experiment 3: Enzyme Inhibition



Serial dilutions of 2,4-PDCA ranging from 3 nM to 300 μM were pre-incubated for 15 min with 0.5 nM JMJD2A. Enzymatic reactions were initiated by the addition of 100 nM biotinylated Histone H3K36me3 peptide substrate plus 2 μM 2OG, 5 μM Fe(II) and 100 μM ascorbate. Enzymatic reactions contain 2% DMSO

Experiment 4: Z'-factor Determination



JMJD2A (0.5 nM) was pre-incubated with or without 100 μM 2,4-PDCA for 15 min. Enzymatic reactions were initiated by the addition of 100 nM biotinylated Histone H3K36me3 peptide substrate plus 2 μM 2OG, 5 μM Fe(II) and 100 μM ascorbate. Enzymatic reactions contain 2% DMSO.



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright © 2010-2011, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

009567_01

Printed in USA

Mar. 2011