


Technical Note

AlphaScreen™

AN004-ASc

ER α Binding Assay



The AlphaScreen™ ER α binding assay has been designed to directly measure the inhibition of ER α binding to estrogen responsive elements (ERE) following exposure to various chemicals. The assay is based on the capture of biotin-ER α and digoxigenin-ERE by streptavidin-donor and anti-digoxin acceptor beads respectively. The AlphaScreen ER α binding assay is specific and reliable. This assay is highly competitive with existing ER α binding assay in terms of ease of use, dynamic range, signal-to-noise ratio (SNR) and time to completion.

The AlphaScreen ER α binding assay is a highly sensitive, homogeneous and non radioactive screening assay. The assay is miniaturized, fully automatable and can be performed in 1 hour.

Technology from:



Packard Instrument Company
800 Research Parkway, Meriden, CT 06450
Tel: 203-639-2598, 1-800-323-1891, Fax: 203-639-2172,
Web site: <http://www.packardinstrument.com>,
E-mail: webmaster@packardinstrument.com

Distributed by:



Example #1:

Competition of digoxigenin-ERE-1 binding to biotin-ER α


developed in Packard 384-well OptiPlate™ microplates (product # 6005214)

Reagents

1. Anti-digoxin-acceptor beads (5 mg/ml in 25 mM HEPES pH 7.4) dilute to 100 μ g/mL in assay buffer
 2. Streptavidin-donor beads (5 mg/mL in 25 mM HEPES pH 7.4) dilute to 100 μ g/mL in assay buffer
 3. Digoxigenin-ERE-1 (50 μ M in PBS) dilute to 7.5 nM in assay buffer with anti-digoxin-acceptor dilution
 4. Biotin-ER α (10 μ M in 25 mM HEPES pH 7.4, 0.1% CHAPS) dilute to 10 nM with anti-digoxin-acceptor beads/ digoxigenin-ERE-1 mix
 5. ERE-1 or ERE-2 analogs (100 μ M in PBS) dilute to 1 nM – 100 μ M in assay buffer
- Assay buffer: PBS + 0.1% BSA

Protocol:

1. Add: 5 μ L ERE-1 or ERE-2
15 μ L digoxigenin-ERE-1/
biotin-ER α -anti-digoxigenin
acceptor beads mix
– incubate 45 minutes at RT
2. Add 5 μ L donor beads
– incubate 15 minutes at RT
3. Read plate



Example #2:

Saturation of ERE-1 by ER α


developed in Packard 384-well OptiPlate™ microplates (product # 6005214)

Reagents:

1. Anti-digoxin-acceptor beads (5 mg/ml in 25 mM HEPES pH 7.4) dilute to 50 μ g/mL in assay buffer
 2. Streptavidin-donor beads (5 mg/mL in 25 mM HEPES pH 7.4) dilute to 100 μ g/mL in assay buffer
 3. Digoxigenin-ERE-1 (50 μ M in PBS) dilute to 7.5 nM in assay buffer with anti-digoxigenin acceptor beads dilution
 4. Biotin-ER α (10 μ M in 25 mM HEPES pH 7.4, 0.1% CHAPS) dilute to 2.5 - 40 nM in assay buffer
- Assay buffer: PBS + 0.1% BSA

Protocol:

1. Add: 5 μ L biotin-ER α , 15 μ L digoxigenin-ERE-1/acceptor beads
– incubate 45 minutes at RT
2. Add 5 μ L donor beads
– incubate 15 minutes at RT
3. Read plate



AlphaScreen™ and OptiPlate™ are trademarks of Packard BioScience Company or its subsidiaries in the United States and/or other countries.