

LANCE *Ultra* EZH2 Histone H3-Lysine 27 N-methyltransferase Assay

U-TRF #40

LANCE® *Ultra*

Authors

Julie Blouin
Mathieu Arcand
Mireille Caron
Anne Labonté
Claire Normand
Lucille Beaudet
Jaime Padrós

PerkinElmer, Inc.
Montreal, QC
Canada, H3J 1R4

This LANCE *Ultra* immunodetection assay measures the dimethylation of a biotinylated Histone H3 (21-44) peptide at lysine 27.

Europium-anti-di/mono-methyl-Histone H3 Lysine 27 (H3K27me₂-1) Antibody

- TRF0406-D: 10 µg, 1,562 assay points*
- TRF0406-M: 100 µg, 15,625 assay points*

*40 fmol/assay point

Peptidic Substrate Sequence:

ATKAARKSAPATGGVKKPHRYRP-GG-K(Biotin)-OH

LANCE *Ultra* Assays

LANCE *Ultra* time-resolved fluorescence resonance energy transfer (TR-FRET) assays use a proprietary europium chelate donor dye, W1024 (Eu), together with *ULight*[™], a small molecular weight acceptor dye with a red-shifted fluorescent emission.

In this technical note, we present the optimization of an epigenetic enzymatic assay using a biotinylated Histone H3-derived peptide as substrate. The modified peptide is captured by the Eu-labeled antibody (Eu-Ab) and *ULight*-Streptavidin (*ULight*-SA), which brings the Eu donor and *ULight* acceptor dye molecules into close proximity. Upon irradiation at 320 or 340 nm, the energy from the Eu donor is transferred to the *ULight* acceptor dye which, in turn, generates light at 665 nm. The intensity of the light emission is proportional to the level of biotinylated substrate modification.




Figure 1. Schematic representation of the LANCE *Ultra* detection of a modified histone peptide.

Development of an EZH2 Histone H3-Lysine 27 N-methyltransferase Assay:

Reagents needed for the assay:


Europium-anti-di/mono-methyl-Histone H3

Lysine 27 (H3K27me2-1) Antibody	PerkinElmer # TRF0406
LANCE Ultra ULight-Streptavidin	PerkinElmer # TRF0102
Histone H3 (21-44) peptide, biotinylated	AnaSpec # 64440
LANCE Detection Buffer, 10X	PerkinElmer # CR97-100
EZH2/EED/SUZ12/RbAp48/AEBP2 Complex	BPS BioScience # 51004
White opaque OptiPlate™-384	PerkinElmer # 6007299
TopSeal™-A films	PerkinElmer # 6005185
Sinefungin	Sigma # S8559
S-(5'-Adenosyl)-L-methionine chloride (SAM)	Sigma # A7007
Poly-L-lysine solution – 0.1% (w/v)	Sigma # P8920

SAM is prepared at 30 mM in 5 mM H₂SO₄/10% ethanol (v/v) in H₂O, aliquoted and stored at -80 °C.


Assay Buffer: 50 mM Tris-HCl pH 9.0, 50 mM NaCl, 1 mM DTT, 0.01% Tween-20 and 0.01% BSA.

Experiment 1: Enzyme Titration and Time-Course



Enzymatic progress curves were performed by incubating EZH2 complex at concentrations ranging from 75 to 200 ng/well with 500 nM biotinylated Histone H3 (21-44) peptide substrate and 30 μM SAM. Reactions were stopped by the addition of poly-L-lysine at indicated times. The Detection Mix was then added and signal was read after 60 min. A 180 min reaction time using 150 ng/well enzyme was selected for all subsequent experiments.

Experiment 2: SAM Titration




Serial dilutions of SAM ranging from 3 nM to 300 μM were added to 150 ng/well EZH2 complex and 500 nM biotinylated Histone H3 (21-44) peptide substrate. A 3 μM SAM concentration was selected for subsequent experiments.

PerkinElmer, Inc.
940 Winter Street
Waltham, MA 02451 USA
P: (800) 762-4000 or
(+1) 203-925-4602
www.perkinelmer.com

Standard Protocol


- Dilute EZH2 enzyme complex, SAM, sinefungin (inhibitor) and biotinylated peptide substrate in Assay Buffer just before use.
- Add to the wells of a white OptiPlate-384:
 - 2.5 μL of enzyme (4X)
 - 2.5 μL of inhibitor (4X) or Assay Buffer
 - 5 μL of biotinylated Histone H3 (21-44) peptide/SAM Mix (2X).
For SAM titration, add SAM dilutions independently of substrate.
- Cover the plate with TopSeal-A film and incubate at room temperature (RT).
- Prepare a 4X Detection Mix by diluting the Eu-Ab to 8 nM and ULight-Streptavidin to 200 nM in 1X LANCE Detection Buffer (final concentrations of 2 nM and 50 nM, respectively, in 20 μL total assay volume).
- Prepare a 4X Stop Solution containing 0.0004% of poly-L-lysine in 1X LANCE Detection Buffer (final concentration of 0.0001% in 20 μL total assay volume).
 - 5 μL of poly-L-lysine Stop Solution and incubate 5 min at RT
 - 5 μL of Detection Mix
- Cover with TopSeal-A film and incubate for 60 min at RT.
- Remove the TopSeal-A film and read signal with the EnVision® Multilabel Reader in TR-FRET mode (excitation at 320 or 340 nm & emission at 665 nm).

Experiment 3: Enzyme Inhibition



Serial dilutions of sinefungin ranging from 300 nM to 30 mM were pre-incubated for 10 min with 150 ng/well EZH2 complex. Enzymatic reactions were initiated by the addition of 500 nM biotinylated Histone H3 (21-44) peptide substrate plus 3 μM SAM. Enzymatic reactions contain 1% DMSO.

Experiment 4: Z'-factor Determination



EZH2 complex (150 ng/well) was pre-incubated with or without 10 mM sinefungin for 10 min. Enzymatic reactions were initiated by the addition of 500 nM biotinylated Histone H3 (21-44) peptide substrate plus 3 μM SAM. Enzymatic reactions contain 1% DMSO.



For a complete listing of our global offices, visit www.perkinelmer.com/ContactUs

Copyright © 2010-2011, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.