

## **GTP $\gamma$ S assay (SPA method) - Sample protocol**

---

[ES-010-M membranes](#) (Adenosine A1 receptor)

*\*Please note that buffer composition, concentration of membrane, concentration of GDP, and concentration of beads will need to be optimized for different receptors.*

### **Reagent preparation**

- Assay buffer 20 mM HEPES pH 7.4; 100 mM NaCl, 10  $\mu$ g/ml saponin, 1 mM MgCl<sub>2</sub>
- Membranes diluted in assay buffer to give 250  $\mu$ g/ml (2.5  $\mu$ g/10  $\mu$ l), keep on ice.
- GDP make a 100  $\mu$ M solution in assay buffer (final concentration of 10  $\mu$ M)
- Beads PVT-WGA (PerkinElmer, RPNQ001), dilute in assay buffer at 50 mg/ml (0.5 mg/10  $\mu$ l)
- GTP $\gamma$  <sup>35</sup>S (PerkinElmer, NEG030X), diluted in assay buffer to give ~25.000 dpm/10 $\mu$ l
- Ligands CPA (Sigma, C-8031), CCPA (Tocris, 1705), GR 79236 (Tocris, 1957), DPCPX (Phoenix, 439) and 1,3-dipropyl-8-phenylxanthine (Tocris, 486); all diluted in assay buffer
- 96 well plate Optiplate™ (PerkinElmer, 6005299)
- Reader TopCount® (PerkinElmer)

### **Make master mixes**

1. Mix the membranes and GDP (volume:volume) and incubate for at least 15 min. on ice
2. Mix the GTP $\gamma$ -<sup>35</sup>S and the beads (volume:volume) just before starting the reaction

### **Add reagents to OptiPlate**

Add sequentially in the order indicated:

3. 50  $\mu$ l of ligand
4. 20  $\mu$ l of the membranes:GDP mix
5. 10  $\mu$ l of buffer\*
6. 20  $\mu$ l of the GTP $\gamma$  <sup>35</sup>S:beads mix

\*For antagonists testing, 10  $\mu$ l of the reference agonist (CPA) are added, at a concentration corresponding to the EC<sub>80</sub> (9.00 nM), instead of 10  $\mu$ l of buffer.

### **Incubation and count**

7. Cover the plate with a TopSeal-A™
8. Shake on an orbital shaker for 2 min
9. Incubate for 1h at RT
10. Centrifuge for 10 min. at 2000 rpm
11. Incubate for 1h at RT
12. Count for 1 min